

The effect of eSight Eyewear on visual function




Walter Wittich, University of Montreal & CRIR
Marie-Céline Lorenzini, University of Montreal & CRIR
Judith E. Goldstein, Johns Hopkins University
Samuel N. Markowitz, University of Toronto
Beatriz E. Patino, University of Toronto
Kristen Lindeman, Johns Hopkins University
Sonya Braudway, Center for Retina and Macular Disease
Scott A. Garter, Lighthouse for the Blind of the Palm Beaches
Lindsay Gotsay, University of Michigan
Ashley Howson, University of Michigan
Michael Tolentino, University of Central Florida
Thiran Jayasundera, University of Michigan
Sophia Reyes, University of Toronto
Gislin Dagnelle, Johns Hopkins University


Disclosure


- This trial is funded by **eSight**
- Marie-Celine Lorenzini holds a peer-reviewed Mitacs Accelerate PhD fellowship, co-funded by Mitacs and eSight



Head-Mounted Displays (HMDs): A History

1990s  Low Vision Enhancement System (LVES)

2000s  Joint Optical Reflective Display (JORDY)


2010s  eSight Eyewear

Things to know about HMDs

- Benefit of hands-free magnification, control over contrast & brightness (Culham, et al., 2004; Geruschat, et al., 1999)
- Initial challenges with size, weight, limited magnification & auto focus (Li et al., 2002)
- Reading can be faster for longer duration and at smaller print sizes than with optical devices (Owsley, et al., 2009; Papageorgiou, et al., Trauzettel-Klosinski, 2010)

eSight Eyewear

- 1) User spots object of interest
- 2) Capture live video
- 3) Digitally optimize the video for the individual in real time
- 4) Project optimized image onto retina using near-to-eye displays




Contrast Enhanced, Magnified

eSight, 2nd generation, 2013

- High resolution 30 fps video camera
- Two SVGA OLED displays (800x600 px)
- Embedded computer similar to iPad and high-end smart phones
- Proprietary image processing circuits and operating system
- Head-borne mass ~200g
- Visual field width 30 deg; diagonal 37.5 deg (3x4 ratio)
- Max magnification 14x

- A full-color digital image is displayed in real-time on high-resolution near-to-eyes organic-LED (OLED) screens




eQUEST Study Objective

- To assess the **impact** of eSight Eyewear on the **functional visual performance** in participants with various low vision conditions.

Method

- Sept 2015 – Dec 2016
- 6 data collection sites










- Johns Hopkins University Wilmer Eye Institute
- University of Michigan Kellogg Eye Center
- The Center for Retina and Macular Disease, FL
- Bascom Palmer Eye Research Institute, FL
- University of Toronto, Toronto Western Hospital
- Université de Montréal, School of Optometry

Participant Characteristics


- Inclusion criteria
 - Age 13 - 75
 - BCVA 20/60 – 20/400
 - Stable central visual loss
 - Min. central field of 20°
 - Ready to wear device in public (e.g. work)
 - Motivated, alert
 - Pass MOCA-Blind cognitive assessment
- N = 51 (30 M, 21 F)
- Mean Age: 48 (SD: 17, Range 13-75)
- Mean Distance VA: 20/178 Range 20/63 - 20/800
- Diagnoses:

– Stargardt's	11
– Leber's	7
– Dry AMD	7
– Optic atrophy	6
– Retinitis Pigmentosa	4
– Rod-cone Dystrophy	3
– Diabetic Retinopathy	1
– Glaucoma	1
– Other (e.g. degen. myopia)	11






Reading Ability



↑ **Endpoint** Significant improvement in reading of **2.76 logits**, $p < .001$
 No site-specific differences

● **Baseline**


Visual Information Items



↑ **Endpoint** Significant improvement in visual information items of **1.08 logits**, $p < .001$
 No site-specific differences

● **Baseline**


Visual Motor Skill Items



↑ **Endpoint** Significant improvement in ADLs of **0.48 logits**, $p < .001$
 No site-specific differences


● **Baseline**

Mobility Items




↑ **Endpoint** No change in mobility items, 0.05 logits, $p = .86$
 No site-specific differences

● **Baseline**



- ### Overall
- eSight Eyewear immediately improved
 - distance acuity
 - reading acuity
 - contrast sensitivity
 - VFQ 48 domains of reading, visual info & visual motor skills
 - We did not observe any significant improvement
 - in reading speed
 - Only reading acuity improved further after 3 months of training & practice
 - Pilot data (not presented today) indicate benefit on
 - face recognition

- ### Meanwhile...
- 3rd-generation eSight device released Feb. 14th 2017
 - Parallel studies:
 - Changes in device use over time
 - Benefit of eSight for pianists while reading sheet music
 - Next steps: Which types of ADLs are most profoundly impacted?
- 

- ### Questions?
- Thank you / Merci:
 - Rob Hilkes – Liaison to eSight
 - Robert Massof – Rasch Analysis
- 

Walter Wittich
 @WalterWittich
 Researcher focusing on the #rehabilitation of individuals with sensory impairment, and on the use and usability of #assistivetechnology
 ▼ Montréal, Québec
opto.umontreal.ca/wittichlab/ind...
walter.wittich@umontreal.ca